

---

**Status:** EU Directives are being published on this site to aid cross referencing from UK legislation. After IP completion day (31 December 2020 11pm) no further amendments will be applied to this version.

---

## [<sup>F1</sup>ANNEX II

### TEST SCHEME FOR DIAGNOSIS, DETECTION AND IDENTIFICATION OF *RALSTONIA SOLANACEARUM* (SMITH) YABUUCHI *ET AL.*

---

#### Textual Amendments

- F1** Substituted by [Commission Directive 2006/63/CE of 14 July 2006 amending Annexes II to VII to Council Directive 98/57/EC on the control of \*Ralstonia solanacearum\* \(Smith\) Yabuuchi et al.](#)

---

*Status: EU Directives are being published on this site to aid cross referencing from UK legislation. After IP completion day (31 December 2020 11pm) no further amendments will be applied to this version.*

---

## Appendix 4

### Buffers for test procedures

#### 1. Buffers for extraction procedure

##### 1.1. Extraction buffer (50 mM phosphate buffer, pH 7,0)

This buffer is used for extraction of the bacterium from plant tissues by homogenisation or shaking.

Na <sub>2</sub> HPO <sub>4</sub> (anhydrous)	4,26 g
KH <sub>2</sub> PO <sub>4</sub>	2,72 g
Distilled water	1.00 l

Dissolve ingredients, check pH and sterilise by autoclaving at 121 °C for 15 min.

Additional components may be useful as follows:

	<b><i>Purpose</i></b>	<b><i>Quantity (per l)</i></b>
Lubrol flakes	Deflocculant <sup>a</sup>	0,5 g
DC silicone antifoam	Anti-foam agent <sup>a</sup>	1,0 ml
Tetrasodium pyrophosphate	Anti-oxidant	1,0 g
Polyvinylpyrrolidone-40000 (PVP-40)	Binding of PCR inhibitors	50 g

<sup>a</sup> For use with homogenisation extraction method]

---