

## SCHEDULE 2

### METHODS OF ANALYSIS

#### 6.

#### *DETERMINATION OF CYANAMIDE NITROGEN*

#### **4.3 REAGENTS**

4.3 Ammoniacal silver solution, according to Tollens, freshly prepared: mix 500 ml silver nitrate solution (10 g per 100 ml) with 500 ml ammonia solution (4.2).

Do not expose unnecessarily to light, heat or air.

*Safety precaution:* when handling ammoniacal silver nitrate solution safety goggles must be worn.