

Status: This is the original version (as it was originally made). This item of legislation is currently only available in its original format.

SCHEDULE 2

METHODS OF ANALYSIS

8b.

DETERMINATION OF DIFFERENT FORMS OF NITROGEN IN THE SAME SAMPLE — IN THE ABSENCE OF CYANAMIDE NITROGEN

Total soluble nitrogen

3.2 *Total soluble nitrogen except nitric nitrogen*, by Kjeldahl digestion after elimination of nitric nitrogen in acid medium by means of ferrous sulphate; ammonia is determined as described in Method 2.