

ANNEX

DETERMINATION OF DICHLOROMETHANE AND 1,1,1-TRICHLOROETHANE

5. APPARATUS

- 5.1. Usual laboratory apparatus.
- 5.2. Gas chromatograph fitted with a thermal conductivity detector.
- 5.3. Transfer bottle, 50 to 100 ml (see sampling method 5.3)⁽¹⁾.
- 5.4. Pressure gas-syringe, 25 or 50 µl (see sampling method 5.4.2.2)⁽¹⁾.

Status: EU Directives are being published on this site to aid cross referencing from UK legislation. After
IP completion day (31 December 2020 11pm) no further amendments will be applied to this version.

- (1) [OJ No L 383, 31. 12. 1980, p. 27.](#)