

## SCHEDULE 2

### Testing Methods

## PART 1

### METHOD FOR THE ISOLATION OF *CLOSTRIDIUM PERFRINGENS*

#### Control Tests

- 16.** Control tests shall be carried out each day that a test is initiated using—
- (a) *Clostridium perfringens* no more than seven days old at the time of use;
  - (b) *Escherichia coli* NCTC 10418(1) or equivalent not more than seven days old at the time of use; and
  - (c) rendered animal protein which is free of *Clostridium perfringens*.